

# SciCan OPTIM<sup>®</sup>

## One-Step Cleaner-Disinfectants



**ON EPA LIST-N OF DISINFECTANTS  
TO USE AGAINST COVID-19**

The global COVID-19 pandemic has affected supply chains in many different industries. In response to the record-breaking demand for OPTIM disinfectants, we have ramped up production. We continue to produce and ship OPTIM Wipes and Solution to dental dealers every week.

The global supply of raw materials used to manufacture pre-saturated wipe products has been substantially constrained as a result of the pandemic. This means you might be experiencing longer than usual lead times on these items.

If you cannot obtain OPTIM Wipes, consider using OPTIM Liquid – **available in both 1 Gallon and 32 oz. bottles**. After all, our Wipes and Liquid use the same formula featured on the Health Canada and EPA Lists of disinfectants recommended against COVID-19. If you cannot obtain OPTIM 1 from your dealer, ask them for OPTIM 33TB.



## USING OPTIM LIQUID

1. Spray SciCan dry wipe or equivalent dry wipe\* with OPTIM solution.
2. If surface is not heavily soiled, clean and disinfect surface with one wipe and allow surface to remain wet for the contact time indicated on the label.
3. If surface is visibly soiled, clean surface with one wipe and discard. Use a second wipe to disinfect and allow surface to remain wet for the contact time indicated on the label.

\*If you cannot obtain SciCan Dry Wipes, you may use paper towel or synthetic towels (e.g. towels made of polyethylene, polyolefin, and polypropylene). Always spray the towel instead of spraying directly onto surfaces.




# OPTIM F.A.Q.



## Is it OK to place disinfectant solution in a container with 4X4 gauze for use on dental equipment?

Use of pre-soaked gauze or paper towel is not recommended. However, it is acceptable to spray OPTIM onto 2x2 or 4x4 gauze or paper towels for immediate use. Gauze is naturally brown in color. The only way to make it white is to heavily bleach it. Presoaking gauze in a disinfectant solution implies putting two competing chemistries head to head in a container which is never a good idea. All gauze tends to absorb too much product and can cause pooling on surfaces, which is inadvisable. If gauze is used for wiping, surfaces should be wiped dry once the contact time has elapsed.

## What is OPTIM's contact time for Novel Coronavirus (COVID-19) and what does that mean?



The contact time is the length of time a disinfectant must be wet on a surface to achieve the desired efficacy result. OPTIM 33TB has a one-minute virucidal contact time and OPTIM 1 has a 30-second virucidal contact time. That means that the surface must remain wet for only 1 minute / 30 seconds respectively in order to inactivate viruses including Novel Coronavirus (COVID-19).

## Can OPTIM be used to clean surfaces?

Yes, OPTIM has been tested as a one-step cleaner disinfectant, meaning it cleans your surfaces as it disinfects. Please refer to the product label or the OPTIM Quick Reference Guide available at [www.scican.com](http://www.scican.com) to learn how to use our product as a one-step cleaner-disinfectant.

## Which surfaces can be cleaned and disinfected by OPTIM?

OPTIM was designed for use on hard, non-porous surfaces. This includes dental chairs, stools, lights and counter tops if used according to our instructions for use on the label. Please refer to the Instructions for Use on the product label.



## Can OPTIM be used to sanitize and disinfect children's toys in dental clinics?

Yes, OPTIM has been classified as non-hazardous according to the GHS requirement. As children's toys do have the possibility to come in contact with children's mouths, the toy must be rinsed with potable water once disinfected.

## OPTIM Ordering Information

### Pre-Saturated Wipes

OPT1-W12	OPTIM 1, Ready-to-Use Wipes, 6" x 7" (12 x 160 sheets)
OPT33-W12	OPTIM 33TB, Ready-to-Use Wipes, 6" x 7" (12 x 160 sheets)

### Liquid

OPT1-4X1	OPTIM 1, Ready-to-Use Liquid (4 x 1G)
OPT1-12X32	OPTIM 1, Ready-to-Use Liquid (12 x 32 oz.)
OPT33-4X1	OPTIM 33TB, Ready-to-Use Liquid (4 x 1G)
OPT33-12X32	OPTIM 33TB, Ready-to-Use Liquid (12 x 32 oz.)

### Accessories

FSA01	Sprayer for 1 L Bottle
SCI33R18	SciCan Dry Wipes, 6" x 7" (18 x 160 sheets)
SCI33R18XL	Large SciCan Dry Wipes, 10" x 10" (18 x 60 sheets)

### Contact Us

#### COLTENE SciCan

701 Technology Dr.  
Canonsburg, PA  
15317 / USA

T +1 724-820-1600

TF +1 800-572-1211

[customerservice@scican.com](mailto:customerservice@scican.com)

[www.scican.com](http://www.scican.com)

 **COLTENE**

OPTIM is a registered trademark of SciCan Ltd.  
OPTIM 33TB is Manufactured by: Virox Technologies Inc. | 2770 Coventry Rd, Oakville, ON L6H 6R1 | Phone: 905-813-0110  
OPTIM 1 is Manufactured for: SciCan Ltd. 1440 Don Mills Road, Toronto, ON, M3B 3P9