

Optim 1 (EU)

Conforms to Regulation (EC) NO. 1907/2006 (REACH), Annex II (453/2010) - Europe

SECTION 1. IDENTIFICATION

Product Identifier Optim 1 (EU)
Recommended Use Disinfectant cleaner
Manufacturer SciCan Ltd, 1440 Don Mills road, Toronto, ON, M3B 3P9, 1-800-667-7733
Emergency Phone No. SciCan Ltd, 1-800-667-7733
SDS No. 053522

SECTION 2. HAZARD IDENTIFICATION

GHS Classification

Not classified under any GHS hazard classes.

GHS Label Elements

Signal Word: None.

Hazard Pictogram: None.

Hazard Statement(s): None.

Other Hazards

This product contains no substances which at their given concentration are considered to be hazardous to health.

SECTION 3. COMPOSITION/INFORMATION ON INGREDIENTS

Mixture:

Chemical Name	CAS No.	%	Other Identifiers
Hydrogen peroxide	7722-84-1	0.5	

Notes

Active ingredient is listed above.

SECTION 4. FIRST-AID MEASURES

First-aid Measures

Inhalation

No specific first aid measures are required.

Skin Contact

No specific first aid measures are required.

Eye Contact

Flush with cool water. Remove contact lenses, if applicable, and continue washing.

Obtain medical attention if irritation develops or persists.

Ingestion

No specific first aid measures are required.

SECTION 5. FIRE-FIGHTING MEASURES

Extinguishing Media

Suitable Extinguishing Media

Not combustible. Use extinguishing agents compatible with product and suitable for surrounding fire.

Unsuitable Extinguishing Media

None known.

Specific Hazards Arising from the Product

None known.

Special Protective Equipment and Precautions for Fire-fighters

As in any fire, wear self-contained breathing apparatus pressure-demand, MSHA/NIOSH (approved or equivalent) and full protective gear.

Product Identifier: Optim 1 (EU)
 SDS No.: 053522
 Date of Preparation: January 22, 2019

SECTION 6. ACCIDENTAL RELEASE MEASURES

Personal Precautions, Protective Equipment, and Emergency Procedures

Use the personal protective equipment recommended in Section 8 of this safety data sheet.

Environmental Precautions

Before attempting clean-up, refer to hazard data. Prevent large spills from entering surface water ways.

Methods and Materials for Containment and Cleaning Up

Small spills or leaks (< 5 gallons) can be flushed into sewage system in accordance with local regulations. Diluted product can be flushed to sanitary sewer. Use water rinse for final cleanup.

SECTION 7. HANDLING AND STORAGE

Precautions for Safe Handling

Use good industrial hygiene practices in handling this material. FOR COMMERCIAL AND INDUSTRIAL USE ONLY.

Conditions for Safe Storage

Store in an area that is out of direct sunlight. Avoid storage at elevated temperatures. KEEP OUT OF REACH OF CHILDREN.

SECTION 8. EXPOSURE CONTROLS/PERSONAL PROTECTION

Control Parameters

Chemical Name	ACGIH TLV®		OSHA PEL		AIHA WEEL	
	TWA	STEL	TWA	Ceiling	8-hr TWA	TWA
Hydrogen peroxide	1 ppm		1 ppm			

Appropriate Engineering Controls

No specific ventilation requirements.

Individual Protection Measures

Eye/Face Protection

Not required if product is used as directed.

Skin Protection

Not required if product is used as directed.

Respiratory Protection

Not required if product is used as directed.

SECTION 9. PHYSICAL AND CHEMICAL PROPERTIES

Basic Physical and Chemical Properties

Appearance	Clear colourless liquid.
Odour	Light almond odour.
Odour Threshold	Not available
pH	2.0 - 2.5
Melting Point/Freezing Point	Not available (melting); Not available (freezing)
Initial Boiling Point/Range	Not available
Flash Point	> 93 °C (200 °F)
Evaporation Rate	Not available
Flammability (solid, gas)	Not applicable
Upper/Lower Flammability or Explosive Limit	Not available (upper); Not available (lower)
Vapour Pressure	Not available
Vapour Density (air = 1)	Not available
Relative Density (water = 1)	1.006 at 20 °C
Solubility	Soluble in water
Partition Coefficient, n-Octanol/Water (Log Kow)	Not available
Auto-ignition Temperature	Not available
Decomposition Temperature	Not available
Viscosity	1.09 centistokes at 20 °C (kinematic)

Product Identifier: Optim 1 (EU)

SDS No.: 053522

Date of Preparation: January 22, 2019

SECTION 10. STABILITY AND REACTIVITY

Reactivity

Not reactive.

Chemical Stability

This product is stable.

Possibility of Hazardous Reactions

Hazardous polymerization will not occur.

Conditions to Avoid

Prolonged exposure to high temperatures.

Incompatible Materials

Do not mix with concentrated bleach products.

Hazardous Decomposition Products

None known.

SECTION 11. TOXICOLOGICAL INFORMATION

Likely Routes of Exposure

Inhalation; skin contact; eye contact; ingestion.

Acute Toxicity

LC50 (Inhalation): > 2.08 mg/L

LD50 (Oral): > 5000 mg/kg

LD50 (Dermal): > 5050 mg/kg

Skin Corrosion/Irritation

Not classified under GHS criteria.

Serious Eye Damage/Irritation

Not classified under GHS criteria.

STOT (Specific Target Organ Toxicity) - Single Exposure

Inhalation

Not classified under GHS criteria.

Skin Absorption

Not classified under GHS criteria.

Ingestion

Not classified under GHS criteria.

STOT (Specific Target Organ Toxicity) - Repeated Exposure

Not classified under GHS criteria.

Respiratory and/or Skin Sensitization

Not classified under GHS criteria.

Carcinogenicity

Not classified under GHS criteria.

Reproductive Toxicity

Development of Offspring

Not classified under GHS criteria.

Sexual Function and Fertility

Not classified under GHS criteria.

Effects on or via Lactation

Not classified under GHS criteria.

Germ Cell Mutagenicity

Not classified under GHS criteria.

Interactive Effects

Not available.

Product Identifier: Optim 1 (EU)

SDS No.: 053522

Date of Preparation: January 22, 2019

Page 03 of 04

SECTION 12. ECOLOGICAL INFORMATION

Ecotoxicity

Not available.

Persistence and Degradability

Not available.

Bioaccumulative Potential

Not available.

Mobility in Soil

Not available.

SECTION 13. DISPOSAL CONSIDERATIONS

Disposal Methods

Review the STORAGE and DISPOSAL instructions on product label prior to disposal.

SECTION 14. TRANSPORT INFORMATION

Not regulated under Canadian TDG regulations. Not regulated under US DOT Regulations.

Special Precautions Not applicable

Other Information RID, ADR (Road, Rail): Not regulated under dangerous goods.

IMO/IMDG (Marine): Not regulated under dangerous goods.

IATA/ICAO (Air): Not regulated under dangerous goods.

SECTION 15. REGULATORY INFORMATION

Safety, Health and Environmental Regulations

Canada

Domestic Substances List (DSL) / Non-Domestic Substances List (NDSL)

All ingredients are listed on the DSL/NDSL.

USA

Toxic Substances Control Act (TSCA) Section 8(b)

All ingredients are on the TSCA Inventory or are exempt from TSCA Inventory requirements under 40 CFR 720.

Europe

Ingredients according to EC Detergent Regulations 648/2004

Anionic Surfactants: Less than 5%

Nonionic Surfactants: Less than 5%

SECTION 16. OTHER INFORMATION

HMIS Rating Health - 0 Flammability - 0 Physical Hazard - 0

SDS Prepared By Virox Technologies Inc.

Phone No. (800) 387-7578

Date of Preparation January 22, 2019

Additional Information For an updated SDS please contact the supplier/manufacturer listed on the first page of this document. Information contained herein was obtained from sources considered technically accurate and reliable. While every effort has been made to ensure full disclosure of product hazards, in some cases data is not available and is so stated. Since conditions of actual product use are beyond control of the supplier, it is assumed that users of this material have been fully trained according to the requirement of all applicable legislation and regulatory instruments. No warranty, expressed or implied, is made and manufacturer/supplier will not be liable for any losses, injuries or consequential damages which may result from the use of or reliance on any information contained in this document. The contents of this document have been prepared in accordance with Regulation (EC) No. 1907/2006, Annex II (453/2010) -Europe.

Product Identifier: Optim 1 (EU)

SDS No.: 053522

Page 04 of 04

Date of Preparation: January 22, 2019