

AQUASTAT™ – Perfect for your office needs.

Every drop. Perfect, steam distilled water.



- **Easy-to-use and energy efficient**

Simply plug it in. No special wiring or plumbing required. Process 3.78 litres (1 gallon) of distilled water in only 4.25 hours. Lightweight design includes conveniently removeable boiling chamber for easy maintenance.

- **Safe, clean water from any source**

Whether your water source is spring fed, municipally piped, soft or hard, the AQUASTAT assures you of fresh, safe, clean, steam-distilled water every time.

- **Protect your equipment and remove water contaminants**

Vaporizes tap water to kill and remove 99% of most contaminants that could damage your equipment.



AQUASTAT

AQUASTAT steam distiller produces perfect, steam distilled water for your autoclave – ideal for STATIM®

Save the environment, stop wasting money on expensive bottled water and avoid damaging equipment from poor quality distilled water.

Weight: 4.5 kg/10 lbs
Dimensions (L x W x H):
 406 mm x 245 mm x 384 mm
 16" x 9 5/8" x 15" 1/8"

Power: 230 VAC, 50Hz, 750W
P/N: W10241S (EU Plug)
 W10243S (Swiss Plug)
 W10244S (UK Plug)