

VISTACLEAR™

Direct Water Filtration, Backflow Prevention
& Operatory Cleaning



POWERFUL PROTECTION AND CONVENIENCE, CENTRALIZED

VistaClear provides dental practices with centrally installed water filtration, backflow prevention and a method for convenient operatory cleaning—eliminating the need for dental bottles altogether! This patented, system saves practices time and money while protecting dental equipment.



Saves Time & Money

VistaClear's on-demand submicron water filtration eliminates wait-time and the need for storage tanks and bottle systems. And its centralized, patented cleaning method makes regular waterline maintenance easier and more efficient.



Protects Equipment


With submicron filtration and no harmful chemicals, VistaClear protects expensive handpieces and turbines from the wear-and-tear typical water supplies can inflict. VistaClear filters remove virtually all microbial contaminants that often reside in city water supplies.



VistaClear Specifications

p/n V1000-28-HP

Operating Temperature Range	7.22-37.78°C (45-100°F)
Operating Pressure Range	30-100 psi (2.07-6.89 bar)
Optimal pH Range	6.5 – 8.5
Maximum pH Range	5.5 – 9.5
Maximum Influent Manganese	0.05 ppm
Maximum Influent Hydrogen Sulfide	0.5 ppm
Maximum Influent Chlorine / Chloramine	5 ppm
Maximum Continuous Flow Rate / Filter	0.2 gpm (0.76 Lpm)
Filter Element Service Life	12 months
Replacement Filter Kit	R9720
Dimensions (W x H x D)	40.64 x 60.96 x 14.92cm
Approximate Ship Weight	28 lbs.



Use of a VistaClear system does not, by itself, eliminate the need for proper waterline maintenance. Due to the nature and complex design of dental delivery units, the periodic use of waterline cleaners and antimicrobials is extremely important for proper dental waterline care. For this reason, Vista Research Group, LLC does not make specific bacteria reduction claims for the output water from dental delivery units.



Upgrade for peak performance.

Do you own an older VistaClear system? You can upgrade your original system to the new design for better performance and lower cost of ownership—for approximately the same price as replacing your filters. Contact your favorite dental dealer for more information.

R9280 VistaClear HP Retrofit System



Serve up to 10 operatories.

VistaClear comes standardly equipped in an eight-operator configuration. It's possible to connect up to 10 operatories to one VistaClear system with an available accessory manifold.

DM125-10C VistaClear 10-Op Distribution Manifold

Distributed by:

SciCan Ltd.

1440 Don Mills Rd., Toronto, ON, Canada, M3B 3P9

Tel: 416-445-1600

Toll free: 1-800-667-7733

customerservice@scican.com

Distributed by:

SciCan Inc.

701 Technology Dr., Canonsburg, PA 15317, USA

Tel: 724-820-1600

Toll Free: 1-800-572-1211

uscustomerservice@scican.com

VistaClear™

is manufactured by:

Vista Research Group

1244 County Road 1475

Ashland, OH 44805 USA

Tel: (419) 281-3927

VistaResearchGroup.com