

VistaPure™

Water Purification System

Meet your steri-center's new best friend

The VistaPure Water Purification System provides two grades of high-quality water for steam sterilizers, dental bottles, ultrasonic cleaners and the final-rinse cycle of compatible instrument washers. The system connects directly to the city water supply and sends high-quality water anywhere you need it in the steri-center for easy, on-demand use.



Two grades of the water practices need

For steam sterilizers:

Point-and-shoot wand makes filling sterilizers with autoclave-ready water a breeze. No more lifting heavy jugs of water!

For bottle systems:

Mounted faucet provides non-corrosive bottle-quality (hyperfiltered) water for bottle systems. No more buying bottled water!



Instant Status Monitoring

The VistaPure Water Purification System features an inline TDS meter that monitors water quality — and makes it easy to know when it's time to replace your filters.



No Electricity Required

The VistaPure Water Purification System operates on water and air pressure, so no electricity is required to supply all the high-quality water your office needs!



Environmentally Friendly

The VistaPure Water Purification System doesn't require chemicals, and its efficient design minimizes the amount of wastewater produced during filtration.

Optional Accessories

Instrument Washer Adapter Kit

S3010

Protect expensive instruments and equipment by using the VistaPure Water Purification System's hyperfiltered water for the final-rinse cycle of instrument washers like the HYDRIM® Instrument Washer.

AutoFill for BRAVO Autoclave Chamber

S3030

Replenish the SciCan BRAVO™ Chamber Autoclave's external refill tank with autoclave-quality water automatically — enabling truly hands-free convenience!

STATCLAVE Adapter Kit

S3015

Send autoclave-quality water directly from the VistaPure System to the SciCan STATCLAVE™ G4 Chamber Autoclave for hands-free automatic filling.

AutoFill for STATIM Autoclave

S3020

Deliver autoclave-quality water to your SciCan STATIM® Autoclave automatically. You'll never have to manually fill your STATIM Autoclave again!




BottleFill Remote Water Dispenser

S3050

Offers the flexibility to choose the optimal placement for bottle filling by mounting securely under an upper cabinet anywhere in the steri-center.



Space Requirements



4.25 Gallon Net Storage Tank S9111

STANDARD

12" Diameter x 18.25" H

2.5 Gallon Net Storage Tank S9110

OPTIONAL

11.125" Diameter x 16" H

9.5 Gallon Net Storage Tank S9113

OPTIONAL

15.5" Diameter x 25" H

Manufactured by:
Crosstex International, Inc.
1733 W. Parkside Lane
Phoenix, AZ 85027 / USA
T +1 (800) 531-2221
VistaResearchGroup.com

Distributed by:
SciCan Ltd.
1440 Don Mills Rd
Toronto, ON M3B 3P9 / Canada
TF +1 (800) 667-7733
T +1 (416) 445-1600
customerservice@scican.com

Distributed by:
SciCan US Office
701 Technology Dr.
Canonsburg, PA 15317 / USA
T +1 724-820-1600
F +1 800-572-1211
uscustomerservice@scican.com

VistaPure is a trademark of Crosstex International, Inc. a Cantel Medical Company.
HYDRIM and STATIM are registered trademarks, and STATCLAVE and BRAVO are trademarks of SciCan Ltd.

For more information, visit www.scican.com

 **COLTENE**